

# ORDER FORM

## DYNAMIC ANKLE FOOT ORTHOSIS (DAFO)

### 1. Company address/Company stamp

\_\_\_\_\_

Phone

\_\_\_\_\_

Date

\_\_\_\_\_

Commission

\_\_\_\_\_

\_\_\_\_\_

**E-Mail** (Confirmation of receipt and tracking number)

\_\_\_\_\_

Responsible staff

\_\_\_\_\_

Signature

### 2. User details

**Has been pre-supplied with a product of ours:**  yes  no **If „yes“** \_\_\_\_\_  
S.-Nr./Date

\_\_\_\_\_

Date of birth

\_\_\_\_\_

User height

\_\_\_\_\_

User weight

**Version:**  a pair  only left  only right

**Activity level:**  1  2  3  4

\_\_\_\_\_

Diagnosis:

\_\_\_\_\_

Additional input:

\_\_\_\_\_

\_\_\_\_\_

### Important Notes:

1. Together with the completed order form we need a rectified and **smooth** plaster positive. Alternatively send a rectified CAD-file to [3d@casimo-ot.de](mailto:3d@casimo-ot.de). A draped check-orthosis to define trim lines and velcro strap positions supports a perfectly fitting device.
2. In case the orthosis includes joints, please set the **pivot point** with the recommended tools from the manufacturer.
3. Please mark the position and the trim lines with a permanent marker the check orthosis or with a **lead pencil** on the plaster.

# ORDER FORM

## DYNAMIC ANKLE FOOT ORTHOSIS (DAFO)

### 3. Execution details

**Please note:** In case there is no selection chosen, our team will follow the steps highlighted in „blue“

**Carbon frame:**  yes  no

**Trimming:**  Casimo-trimming  according to your drawings

**Padding:**  malleolus  talus (medial)  complete foot  partial foot

**Padding thickness:**  4  5  6

**Padding color:**  black  blue  yellow  red  white  pink

**Plastic material:**  Copolymer  PP  PE

**Thickness:**  2  3  4

**Plastic color:** \_\_\_\_\_  
(transfer paper, please see [casimo-ot.de](http://casimo-ot.de))

**Carbon color:** \_\_\_\_\_  
(decorative fabric, please see [casimo-ot.de](http://casimo-ot.de))

**Velcro strap:**  elastic  rigid

**Velcro strap color:**  black  white  blue  yellow  red